



QuikGlaze™ - QuikGlaze™ XL & QuikGlaze™ Temp White

Toll Free: 877-64-PRODX
Fax: 877-408-0114
www.alldentalprodx.com

Sealant and Bonding Agent for ALL Temporaries

PRODUCT DESCRIPTION



Dual-purpose **QuikGlaze™** and **QuikGlaze™ Temp White** are light-cured sealants and bonding agents for use with ANY and ALL provisional materials. Fill microscopic voids. Temp will be abrasion and stain resistant. Adds natural looking esthetics to any provisional, regardless of material used. Eliminate the need for polishing.

“Finally a quick, easy, one-step process to simplify and speed up chairside temp polishing. I appreciate a product like **QuikGlaze™** because it allows me to focus on the task at hand...providing my patients with the best possible service!”

Steven Brattesani, DDS

FEATURES

Produces a smooth, esthetic finish

Both a sealant and a bonding agent

Works with *any and all* temporary materials

Brush on, Light cure

QuikGlaze™ Temp White has a built in whitener and brightener

BENEFITS

Eliminates polishing time
Gives provisional restorations a natural looking luster
Extends longevity: stain and abrasion resistant

Multi-purpose product. Finish temps and quickly repair flawed temps: voids or inadequate margins.

Versatility

Ease and speed of use reduces chair time.

Inventory fewer shades of temporary materials
Mix with original **QuikGlaze™** for optimum esthetics

PRODUCT ORDERING INFORMATION

200130

200130-XL (not shown)

200150

QuikGlaze™ 5ml Bottle / Accessories

QuikGlaze™ XL 30ml Bottle / Accessories

QuikGlaze™ Temp White 5ml Bottle / Accessories

Call PRODX toll-free to place your order or call your dental dealer directly.

Michael Miyasaki DDS, LVIM

Dr. Michael Miyasaki is LVI's Vice President of International Operations. A 1987 graduate of USC School of Dentistry, he developed a highly successful reconstruction practice in Sacramento, CA. Following his passion to teach and mentor other dentists, he became associated with LVI in 1996 where he now works full time. Michael practices in the LVI faculty practice, lectures and publishes articles on the latest aesthetic, occlusion and materials available.



All Dental PRODX QuikGlaze


*by Michael Miyasaki DDS, LVIM
Article from LVI VISIONS, Jan/Feb/Mar 2006*



How temporaries often look after we have cleaned them up. They appear rough.



Painting on QuikGlaze



Light curing the QuikGlaze material



Right central with QuikGlaze

Imagine, it's the end of that veneer appointment. You're tired, your patient is tired and you still have to make the temps shine. One option is polishing the temps with the various cups, disks and points in your composite armamentarium. Or another more appealing option is placing a light cured glaze over the top of temporaries and light curing leaving you with a shiny smooth surface free of the sticky oxygen inhibited layer. The temporaries not only look amazingly smooth, but are more stain resistant. This process is fast and easy. Try it as that final step to provide your patients with that great 'wet' look in just a minute's worth of work.

The list of ingredients lists 'multifunctional acrylates' which enables you to use this material as a bonding agent when repairing or adding to provisional materials. Just roughen the surface to be repaired, apply a layer of QuikGlaze, then cure the material and place the repair material on top. After curing you have bonded the layers together.

So this versatile glaze material is a light-cured methyl methacrylate varnish that can be used in the following situations:

- Repair of provisional resin materials
- Seals surfaces of temporary resin materials
- Eliminates time consuming polishing of resin temporaries

In this issue I've discussed... products that are proved performers at LVI. With the hundreds of doctors we train each year these products have proven their worth in the clinic... a fast glazing agent that will save you valuable minutes while providing your patients with knock-out smiles while in their temporaries.